

ApoCure®

- hydrogel dressing



The Hydrogel contains only 30% water, and 70 % absorbent polymer. The ApoCure Hydrogel is permeable to water vapour. The ApoCure Hydrogel creates a moist wound environment and will not adhere to the under laying tissue at dressing changes.

ApoCure is transparent and secures that the wound can be inspected. ApoCure Hydrogel can remain on the treated area up to seven days.

ApoCure Hydrogel allows the exudates to be absorbed into the dressing through vertical uptake. Based on the absorption features of the dressing, - allows the exudate to be bound in the dressing, and the dressing can by that be used under compression therapy for venous leg ulcer treatment.

The Hydrogel coating secure that the dressing doesn't adhere to the wound bed, and by that minimizes the trauma and pain during dressing changes. Studies shows that bacteria will be absorbed by the Hydrogel and by osmotic action destroyed.

Studies shows that patients report a higher degree of pain relieving effect.



Item	Beskrivelse	Size	Dressings/box	Box / Carton
021525	ApoCure Oval with border	4,5 x 7,5 cm	50	8

For more information:

Ph: +45 32 97 15 25 | www.apodan.dk | info@apodan.dk